


POWER SCOPE™ 2



CLASS II CORRECTION SIMPLIFIED

PowerScope 2 delivers simple, efficient Class II correction like you've never seen before. The latest generation of PowerScope features a stronger wire attachment nut and visual activation indicators.

- Fixed one-piece, one-size-fits-all design
- Easy chairside wire-to-wire installation
- No assembly or lab setup required
- No headgear or special band assemblies
- Use with banded or bondable molar tubes
- World-class clinical results



NEW **MAGNETIC FEATURE**

Magnetic sleeve slides onto driver for secure hold of appliance during installation.



Scan for more information.

FEATURES AND BENEFITS

- One-size-fits-all Class II corrector for simple treatment application and easy inventory management
- One-piece, no assembly required, and no need for lab setup
- Internal NiTi spring mechanism delivers 260 grams of force during treatment
- Connects wire-to-wire for quick and easy chairside installation
- Can be used with both bondable or banded molar tubes
- No molar headgear tube required
- Telescoping mechanism will not disengage during treatment, helping avoid emergency visits
- Unique ball and socket joint for unrestricted movement and patient comfort
- Fatigue tested to ensure durability during treatment
- Can be used for unilateral or bilateral correction
- Works with fully bonded cases
- Crimpable shims available for activation of spring during treatment

POWERSCOPE KITS

- 852-016537** 1 Patient Kit
(Includes: 1 Right/1 Left with nuts – *does not include Shims, Driver, or Magnets*)
- 852-016538** 5 Patient Kit
(Includes: 5 Right/5 Left with nuts, 10 each 1 mm, 2 mm & 3 mm Crimpable Shims, 1 Hex Head Driver, and 5 Driver Magnets)

POWERSCOPE ACCESSORIES

- 852-016539** PowerScope 2 — Right Arm (18 mm with nuts)
- 852-016540** PowerScope 2 — Left Arm (18 mm with nuts)
- 852-016541** Replacement Attachment Nuts — 10 pack (5 right and 5 left)
- 855-733** Crimpable Shims 1 mm — 10 pack
- 855-728** Crimpable Shims 2 mm — 10 pack
- 855-738** Crimpable Shims 3 mm — 10 pack
- 855-515** Hex Head Driver
- 852-026533** Driver Magnets — 5 pack
- 855-515** Hex Head Driver



3524 Washington Avenue
Sheboygan, Wisconsin, 53081 USA

Phone: +1 920 457 5051
Email: cs@americanortho.com
Web: americanortho.com



©2017 American Orthodontics Corporation
All Rights Reserved

FOLLOW US ON

